




HT reactor plate, complete, for Professional Reactor

6.2845.110

Reactor plate with reactor for PCR applications at high temperatures ($> 120\text{ }^{\circ}\text{C}$). With PEEK reactor capillaries.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.2845.110

Qt.	Order no.	Description
-----	-----------	-------------

3 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description
6.2057.240	Holder for reactor 62845xxx Holder for fastening a reactor in the column thermostat.

