



Exchange Unit 5 mL

6.3026.150

Exchange Unit with integrated data chip with 5 mL glass cylinder and light protection. PCTFE /PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.


Scope of delivery 6.3026.150

Qt.	Order no.	Description
1 PCS	6.1543.060	Tip / M6 thread

Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



Optional accessories

Order no.	Description	
6.1563.040	Thermostat jacket for 806 Exchange Units For 6.3026.XXX Exchange Units	
6.1576.150	Exchange unit base plate / 5 mL	
6.1608.030	Round glass bottle / 1000 mL / GL 45	
6.1608.040	Bottle PE / 1000 mL / GL 45 For Exchange Units. Bottle for auxiliary solutions.	