



Dosing unit PCR 5 mL

6.3033.150

5 mL dosing unit specially suitable for conveying solutions for post-column derivatization (PCR) in IC.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.3033.150

Qt.	Order no.	Description
1 PCS	6.2739.000	Wrench For tightening connectors



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Optional accessories

Order no.	Description
6.2057.250	Dosino holder for PCR
	Dosino holder for IC instruments for PCR Dosino.

