



Differential amplifier plug F / 230 V EU

6.5104.030

Differential amplifier with power supply.


Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.5104.030

Qt.	Order no.	Description
1 PCS	6.2117.000	Power adapter EU 230 V / 9 V AC With EU mains plug. Input: 220...240 V AC. Output: 9V AC/0.2 A. Used with 587 Conductometer, 588 pH Meter, 627 O ₂ Meter



Differential amplifier with plug F, for Titrados (power supply unit must be ordered separately: 6.2117.000 power unit EU 230 V / 9 V AC)

Also available as a complete set: 6.5104.030 differential amplifier / plug F / 230 V EU

