



IC equipment: Low Volume Inline Dialysis

6.5330.200

Accessory set for rapid Inline Dialysis. For use with the 858 Professional Sample Processor and an additional two-channel peristaltic pump.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.5330.200

Qt.	Order no.	Description
2 PCS	6.2744.070	Pressure screw short
		Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.1608.070

Eluent bottle / 2 L / GL 45



1 PCS

6.1831.160

PEEK capillary, 0.5 mm ID / 70 cm



2 PCS

6.2744.034

Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.




4 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



1 PCS

6.1602.150

Bottle attachment / GL 45 - 3 x UNF 10/32

For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.



8 PCS

6.1803.040

PTFE capillary 0.5 mm i.d. / 1 m

Capillary for sample handling in IC.



2 PCS

6.1826.340

Pump tubing LFL (black/black), 3 stoppers

For the sample solution in inline dialysis.



1 PCS

6.2755.000

Tubing cartridge

For the peristaltic pumps of IC instruments. Contact pressure adjustable.



1 PCS

6.1831.120

PEEK capillary 0.25 mm i.d., 45 cm



1 PCS

6.2744.000

Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



2 PCS

6.1826.320

Pump tubing LFL (orange/yellow) , 3-stopper

For suppressor solutions and for acceptor solutions for inline dialysis and for inline ultrafiltration.



1 PCS

6.2831.010

Plastic tweezers

For dialysis and filtration membranes




1 PCS

6.2057.130

Dialysis cell holder

Holder for the dialysis cell



1 PCS

6.2714.010

Dialysis membrane (cellulose acetate)

Nominal pore size : 0.2 µm. Set of 50 pieces. For inline dialysis



1 PCS

6.2729.200

Low Volume dialysis cell

Including four 6.2744.060 stoppers. For rapid Inline Dialysis.

