



# IQ/OQ Kit for IC

6.5333.000

The IQ/OQ Kit contains all parts and standard solutions required for IQ/OQ in ion chromatography.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 6.5333.000

| Qt.   | Order no.  | Description                                                                                     |
|-------|------------|-------------------------------------------------------------------------------------------------|
| 1 PCS | 6.1803.050 | PTFE capillary 0.5 mm i.d. / 20 cm                                                              |
|       |            | Capillary for inline dialysis and filtration, for Dialysis Unit, IC Filtration Sample Processor |



1 PCS

6.1831.170

OQ capillary

PEEK capillary used in OQ procedures for Professional IC equipment



---

1 PCS

6.2322.030

OQ Standard Solutions for IC

Part of the IQ/OQ Kit for IC (6.5333.000). Cannot be ordered individually.



---

1 PCS

6.2744.040

2 x UNF 10/32 coupling

For connecting 1/16 in. capillaries. For IC instruments



---

1 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



1 PCS

6.2744.070

Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries

