



Grounding equipment for Professional VA systems with automatic addition of solution

6.5334.020

Accessory set for two auxiliary solutions for expanding semiautomated Professional VA systems. Contains one ground adapter for two Dosino instruments and one connection cable.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 6.5334.020

Qt.	Order no.	Description
-----	-----------	-------------

2 PCS

6.2167.000

Ground adapter for Dosino

Connector with M6 thread outer for Dosino and M6 inner for 6.1824.000 tubing set



1 PCS

6.2167.010

Ground cable, 1 m, USB A / plug 2 mm

For the connection of one 6.2167.000 ground adapter to a USB connector.



1 PCS

6.2167.020

Connecting cable, 0.5 m, 2 x plug 2 mm

