




# IC equipment: Waste Collector

6.5336.000

Waste Collector for up to five waste capillaries.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 6.5336.000

| Qt.   | Order no.   | Description                   |
|-------|-------------|-------------------------------|
| 1 PCS | 6.01816.010 | Silicone tubing, 8/12 mm, 2 m |
|       |             | Drainage tubing               |



1 PCS

6.1602.210

Attachment 5 x UNF 10/32 for «Waste Collector»

Attachment for holding up to 5 waste capillaries.



---

1 PCS

6.2057.220

Holder to Waste Collector

Holder for fastening the Waste Collector to the IC instrument.



---

1 PCS

6.2765.000

Container with GL 45 drainage to the "Waste Collector"

Container with drainage to the "Waste Collector".

