



Equipment with 2x 807 Dosing Units for VA and CVS measurements

6.5339.500

Accessory set for 2 auxiliary solutions for expanding VA or CVS systems. Contains 2x 807 Dosing Units 2 mL, 2x 100 mL clear glass bottles, bottle holder for 2 bottles, and connector material. Without Dosinos.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5339.500

Qt.	Order no.	Description
-----	-----------	-------------

2 PCS

6.1608.050

Drying bottle / 100 mL / GL 45



2 PCS

6.1618.060

Thread adapter GL 45 to GL 45

Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange



2 PCS

6.1805.530

FEP tubing / M6 / 2 m

With light and kink protection



2 PCS

6.1819.010

PTFE tube

As gas inlet in washing bottle of VA Stand



1 PCS

6.2055.110

Bottle holder or Dosino holder for two bottles
/Dosinos

Holder for two 250 mL bottles or two Dosinos



2 PCS

6.3032.120

807 Dosing Unit 2 mL

807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



Optional accessories

Order no.	Description
2.800.0010	800 Dosino The 800 Dosino is a drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm). The image shows the 800 Dosino drive unit, which is a white, L-shaped device. A coiled white cable is attached to the unit. Next to the unit is a small white probe with a metal tip.
