



Measuring stand (scTRACE Gold) replacement for the 946 Portable VA Analyzer

6.5340.010

Measuring stand for the scTRACE Gold, complete, comprising a base plate, transparent ring, and measuring head.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.01204.000

Stirrer tip


Stirrer tip for 946 Portable VA Analyzer.



6.01408.000

O-ring 46 x 2.5 mm

For the 946 Portable VA Analyzer measuring stand. Set of 2 pieces.



6.02621.000

4 mm wrench



6.1258.000

scTRACE Gold (set with 4 pieces)

Gold micro-wire electrode on thick film basis for the determination of arsenic in waters. The scTRACE Gold is a combined voltammetric sensor for voltammetry. Gold micro-wire working electrode, Ag/AgCl reference electrode and carbon auxiliary electrode are located on a sensor. Simple in its handling and practically maintenance-free, the scTRACE Gold can for example be used to determine arsenic in water in the trace range. Together with the 6.1241.080 electrode shaft, it can be used with any Metrohm voltammetry measuring stand.

