



Equipment for determination of temperature correction with Biodiesel Rancimats.

6.5616.110

Set for exact temperature adjustment

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.5616.110

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1111.020

Pt100 temperature sensor for stability measuring devices (30 cm)

Pt100 temperature sensor (length 30.0 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with instruments for stability measurement (873 Biodiesel Rancimat and 893 Professional Biodiesel Rancimat).



2 PCS

6.2042.040

Clamp for temperature sensor

Holding clamp for Pt100 temperature sensor for stability measuring instruments



1 PCS

6.2326.000

Silicone oil for stability measuring instruments (50 mL)

Silicone oil for the determination of temperature correction in stability measuring instruments

