




Measuring equipment for checking the temperature

6.5618.000

Measuring equipment for checking the temperature of 774 Oven Sample Processor, 832, 860 Thermoprep.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5618.000

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.2621.120

Hexagon key 1.5 mm



1 PCS

6.1253.010

Measuring insert

For checking the oven temperature of the 774 Oven Sample Processor

