



Transfer equipment, 0.5 - 2 mL, for 789/815

6.5619.000

Complete equipment for the transfer of solutions from 0.5 to 2.0 mL on the Robotic Sample Processor.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.5619.000

Qt.	Order no.	Description
1 PCS	6.1446.180	Stopper, PTFE, SGJ 29, for transfer set
		Stopper for titration and volumetric KF titration in the external titration vessel



1 PCS

6.1448.000

Septum stopper B-14/15



2 PCS

6.1448.020

Septum 16 mm, 5 pieces

Set of 5 pcs.



1 PCS

6.1462.090

Robotic arm with Luer-Lock adapter for 786 Swing
Head, right swinging

Transfer robotic arm for liquid handling



1 PCS

6.1464.320

KF titration vessel / 80-250 mL / coulometric



1 PCS

6.1562.110

Pipetting tubing / 2 mL

Pipetting tubing for liquid transfer



1 PCS

6.1608.090

Bottle / 250 mL / GL 45




1 PCS

6.2001.100

Base plate, transfer set

Base plate for transfer set for liquids to be used with Robotic Sample Processors XL



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2701.040

GL 18 screw cap with hole

Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.



1 PCS

6.2751.120

Safety shield Transfer set to 789/815

Safety shield for Sample Processors with 786 Swing Head and transfer set.



1 PCS

6.2816.100

Needle with vent

With luer lock, for transfer equipment 6.5619.000.



1 PCS

6.2055.110

Bottle holder or Dosino holder for two bottles /Dosinos

Holder for two 250 mL bottles or two Dosinos

