




# Consumable Kit Titration

6.5702.010

FEP tubing connection M6 (13, 25 and 40 cm) and antidiffusion tip for 6.3026.xxx exchange units and 6.3032.xxx dosing units.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 6.5702.010

| Qt.   | Order no.  | Description                             |
|-------|------------|-----------------------------------------|
| 1 PCS | 6.1543.200 | Titration tip / M6 thread               |
|       |            | Titration tip with antidiffusion valve. |



2 PCS                      6.1805.010                      FEP tubing / M6 / 13 cm  
With light and kink protection.



---

2 PCS                      6.1805.080                      FEP tubing / M6 / 25 cm  
With light and kink protection.



---

2 PCS                      6.1805.100                      FEP tubing / M6 / 40 cm  
With light and kink protection



---

1 PCS                      6.1829.010                      FEP aspiration tube M6, 0.25 m  
With M 6 thread. For Dosing Unit



4 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,  
6.1543.150, 6.1543.160, 6.1543.200

