




Kit for liquid samples for MMS 5000, 112 positions

6.7303.000

Kit for liquid samples for MMS 5000 with rack for sample vessels, 100 μ L syringe with replaceable needle and 100 snap-ring bottles.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.7303.300

CIC Liquid Sampler solvent reservoir

Reservoir for rinsing the sample needle on the Autosampler MMS 5000 with liquids kit.

