



Glass vial, 1.5 mL, 100 pieces

6.7303.100

For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.


Optional accessories

Order no.	Description
6.2041.200	Sample rack, 48 x 11 mL, (Al) Sample rack for 889 IC Sample Center, suitable for 48 sample vessels. The following sample vessels can be used: 6.2743.100, 6.2743.107, 6.2743.110, 6.2743.117. Material: Aluminum (Al)



6.7303.110 Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers:
Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889
Sample Center.



6.7303.310 Sample rack 112 positions for liquid samples

For liquid samples.

