



Titration head for CSB

6.9914.137

Titration head for use in the CSB System (MATi 12). The titration head offers three retracted tubings with antidiffusion valve and one for aspiration of the titrated solution. The electrode can be fastened in place in an M10 screw connector. As before, three openings are provided for spray nozzles.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.9914.137

Qt.	Order no.	Description
-----	-----------	-------------

4 PCS

6.1808.000

Connection bushing M6



3 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200

