



Overflow buret (100 mL) for Sample Processor

6.9958.007

Overflow buret for rapid and automatic measuring of 100 mL of sample and for sample transfer into an external titration vessel.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.9958.007

Qt.	Order no.	Description
2 PCS	6.1808.050	Tubing adapter olive / M8 outer
		1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



2 PCS

6.1801.000

PVC tubing / 4 mm / 7 mm / 1 m



1 PCS

6.1819.040

FEP tube

For M8 threaded holes. For 6.1563.010 thermostat jacket of 6.3011.XXX...
6.3026.XXX Exchange Units

