



843 Pump Station (membrane)

2.843.0020

L'843 Pump Station (membrana) possiede due pompe a membrana integrate. Queste possono essere comandate a scelta tramite l'interfaccia direttamente tramite segnali remoti o manualmente con la pressione di un tasto.

Parti incluse 2.843.0020

Qt.	Order no.	Descrizione
1 PCS	1.843.0020	843 Pump Station (membrane)
		L'843 Pump Station (membrana) possiede due pompe a membrana integrate. Queste possono essere comandate a scelta tramite l'interfaccia direttamente tramite segnali remoti o manualmente con la pressione di un tasto.



4 PCS

6.1801.120

Tubo in PVC / 4 mm / 6 mm / 2 m



2 PCS

6.1812.000

Tubo flessibile PTFE 4 mm / 6 mm / 4 m



2 PCS

6.1828.000

Connessione in PVDF

Per recipiente 6.1621.000.



Accessori opzionali

Order no.	Descrizione	
6.2141.260	Remote cable Compact Sample Changer - Dosimat plus and Titrino plus Cable for remote connection between Compact Sample Changer and Dosimat plus and further devices (e.g. Titrino plus or Pump Station).	
6.2141.300	Cavo remoto 843-Sample Processor Cavo remoto per la connessione diretta di una 843 Pump Station alle connessioni della pompa del Sample Processors	