



Balance for 864 Robotic Balance Sample Processor, 115V

2.135.0205


納品範囲 2.135.0205

Qt.	Order no.	説明
1 PCS	1.135.0205	Balance for 864 Robotic Balance Sample Processor, 115V

The balance (Precisa) has a special construction for usage below the sample rack 6.2068.010 of the 864 Robotic Balance Sample Processor. The power supply is for 115V.



オプションアクセサリ

Order no.	説明	
6.2067.100	Scale pan for sample rack 6.2068.010	
6.2068.010	サンプルラック、20 x 120 mL、(ポリプロピレン製) 20 個のサンプルビーカーに適した864 Balance Sample Processor のためのサンプルラック。以下のサンプルビーカーの使用が可能です: 6.1459.300。 合成樹脂: ポリプロピレン (PP)	
6.2711.080	Drip pan for 864 Robotic Balance Sample Processor	