



# Balance for 864 Robotic Balance Sample Processor, 230V

2.135.0206


納品範囲 2.135.0206

| Qt.   | Order no.  | 説明                                                     |
|-------|------------|--------------------------------------------------------|
| 1 PCS | 1.135.0206 | Balance for 864 Robotic Balance Sample Processor, 230V |

The balance (Precisa) has a special construction for usage below the sample rack 6.2068.010 of the 864 Robotic Balance Sample Processor. The power supply is for 230V.



## オプションアクセサリ

| Order no.  | 説明                                                                                                                                                         |                                                                                       |
|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.2067.100 | Scale pan for sample rack 6.2068.010                                                                                                                       |    |
| 6.2068.010 | サンプルラック、20 x 120 mL、(ポリプロピレン製)<br><br>20 個のサンプルビーカーに適した864 Balance Sample Processor のためのサンプルラック。以下のサンプルビーカーの使用が可能です: 6.1459.300。<br><br>合成樹脂: ポリプロピレン (PP) |    |
| 6.2711.080 | Drip pan for 864 Robotic Balance Sample Processor                                                                                                          |  |