



# 705 UV Digester (220 V, 60 Hz)

2.705.0016

納品範囲 2.705.0016

| Qt.   | Order no.  | 説明                             |
|-------|------------|--------------------------------|
| 1 PCS | 1.705.0016 | 705 UV Digester (220 V, 60 Hz) |



1 PCS

1.705.0020

705 UV Digester wet end



---

1 PCS

6.2804.030

UV Mercury vapor lamp



---

12 PCS

6.1446.100

PTFE stopper for 12 mL UV quartz sample vessel



---

1 PCS

6.1446.110

Thermometer stopper for UV quartz sample  
tube 12 mL



1 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



---

1 PCS

6.2041.040

Sample holder for 12 quartz digestion vessels,  
12 mL



---

1 PCS

6.2401.010

Mercury thermometer



---

12 PCS

6.2414.000


Quartz sample vessel 12 mL



1 PCS

6.2745.000

UV protection plate



## オプションアクセサリ

| Order no.  | 説明                                                  |                                                                                       |
|------------|-----------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.1446.150 | Thermometer stopper for UV quartz sample tube 5 mL  |    |
| 6.2041.050 | Sample holder for 16 quartz digestion vessels, 5 mL |    |
| 6.1446.140 | Stopper for 5 mL UV quartz sample vessel            |  |
| 6.2414.010 | Quartz Sample Vessel 5 mL                           |  |