



705 UV Digester (220 V, 50 Hz)

2.705.0017

納品範囲 2.705.0017

Qt.	Order no.	説明
1 PCS	1.705.0017	705 UV Digester (220 V, 50 Hz)



1 PCS

1.705.0020

705 UV Digester wet end



1 PCS

6.2804.030

UV Mercury vapor lamp



12 PCS

6.1446.100

PTFE stopper for 12 mL UV quartz sample vessel



1 PCS

6.1446.110

Thermometer stopper for UV quartz sample tube 12 mL



1 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.2041.040

Sample holder for 12 quartz digestion vessels,
12 mL



1 PCS

6.2401.010

Mercury thermometer



12 PCS

6.2414.000


Quartz sample vessel 12 mL



1 PCS

6.2745.000

UV protection plate



オプションアクセサリー

Order no.	説明	
6.1446.150	Thermometer stopper for UV quartz sample tube 5 mL	
6.2041.050	Sample holder for 16 quartz digestion vessels, 5 mL	
6.1446.140	Stopper for 5 mL UV quartz sample vessel	
6.2414.010	Quartz Sample Vessel 5 mL	