

733 IC Separation Center

2.733.0020

納品範囲 2.733.0020

Qt.	Order no.	説明	
1 PCS	6.1816.000	Silicone tubing 6 mm, 0.5 m	
1 PCS	6.1831.010	PEEK capillary 0.25 mm i.d. / 3 m	
2 PCS	6.2027.030	Column holder 8.5 mm	

2 PCS

6.2027.040

Column holder 11.3 mm



4 PCS

6.2620.060

Coupling



1 PCS

6.2621.000

Adjustable wrench



1 PCS

6.2621.050

1/4 in. wrench



1 PCS

6.2621.060

5/16 in. wrench



2 PCS

6.2744.010

Pressure screw 5x



2 PCS

6.2744.020

Luer/UNF coupling




2 PCS

6.2816.020

Syringe 10 mL with Luer connection



オプションアクセサリ

Order no.	説明	
6.2620.110	Sampling loop 20 μ L	
6.2620.120	Sampling loop 100 μ L	
6.1825.210	PEEK sample loop 20 μ L	
6.1825.230	PEEK sample loop 10 μ L	



6.5904.030 ローター、6/2 インジェクタ

IC システムの 6/2 インジェクタ用ローター。



6.5904.050 内部ループを伴うインジェクタ 0.25 µL

内部ループを伴うフラットプレートインジェクタ 0.25 µL



6.5904.060 ローター（内部ループ、0.25 µL）

4ポートフラットプレートインジェクタのためのローター
(6.5904.050)、内部ループ 0.25 µL



6.5904.070

ローター (内部ループ、1.0 μ L)

4ポートフラットプレートインジェクタのためのローター
(6.9956.001)、内部ループ 1.0 μ L

