

Pump Unit 772 "Aspiration"

2.772.0020

納品範囲 2.772.0020

Qt.	Order no.	説明	
1 PCS	6.1446.160	Thermometer stopper, B-14/15, oblique insertion	
1 PCS	6.1543.170	Aspiration tip / M8 thread	
1 PCS	6.1621.000	Container 10 L	

1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1820.050

M6 / Olive adapter



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m



1 PCS

6.1828.000

PVDF connection nipple



1 PCS

6.2621.140

Hexagon key 2.5 mm



オプションアクセサリー

Order no.	説明	
6.1826.160	Viton pump tubing	
