



# 772 Pump Unit – aspirate

2.772.0120

納品範囲 2.772.0120

| Qt.   | Order no.  | 説明            |
|-------|------------|---------------|
| 1 PCS | 1.772.0110 | 772 Pump Unit |

Peristaltic pump for direct connection to Sample Processors.



1 PCS

6.1446.160

Thermometer stopper, B-14/15, oblique  
insertion



---

1 PCS

6.1543.170

Aspiration tip / M8 thread



---

1 PCS

6.1621.000

Container 10 L



---

2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm



---

1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



---

2 PCS

6.1820.050

M6 / Olive adapter



---

1 PCS

6.1820.070

Coupling 4 mm / 6 mm / M8 inner



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m



---

1 PCS

6.1828.000

PVDF connection nipple



オプションアクセサリ

| Order no. | 説明 |
|-----------|----|
|-----------|----|

6.1826.160      Viton pump tubing



6.2160.010      Adapter cable / 2 m / M8 socket / 2 x plug B

