



772 Pump Unit – rinse

2.772.0130

納品範囲 2.772.0130

Qt.	Order no.	説明
1 PCS	1.772.0110	772 Pump Unit

Peristaltic pump for direct connection to Sample Processors.



1 PCS

6.1621.000

Container 10 L



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



3 PCS

6.1805.060

FEP tubing / M6 / 60 cm



1 PCS

6.1808.170

Distributor for Sample Processors for rinsing and aspirating



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1820.050

M6 / Olive adapter



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m



1 PCS

6.1828.000

PVDF connection nipple



3 PCS

6.2740.020

Spray nozzle



オプションアクセサリ

Order no.	説明
-----------	----

6.1826.160 Viton pump tubing



6.2160.010 Adapter cable / 2 m / M8 socket / 2 x plug B

