



# 786 Swing Head with transfer head left

2.786.0010

納品範囲 2.786.0010

| Qt.   | Order no.  | 説明                                     |
|-------|------------|----------------------------------------|
| 1 PCS | 1.786.0010 | 786 Swing Head with transfer head left |



1 PCS

6.1462.030

Robotic arm with transfer head for 786 Swing  
Head, left swinging



1 PCS

6.1562.100

Pipetting tubing / 10 mL



1 PCS

6.1823.010

Guide shaft to pipetting tubing



1 PCS

6.2621.120

Hexagon key 1.5 mm

