



786 Swing Head with titration swing arm, left- or right-swinging

2.786.0030

納品範囲 2.786.0030

Qt.	Order no.	説明
1 PCS	1.786.0010	786 Swing Head with transfer head left



1 PCS

6.1462.050

Robotic arm with titration head for 786 Swing
Head, left or right swinging



1 PCS

6.1808.040

Thread adapter M6 outer / M8 inner



1 PCS

6.1909.030

Stirring propeller / 95 mm



2 PCS

6.2042.030

Buret tip clip (rod stirrer)



1 PCS

6.2621.120

Hexagon key 1.5 mm

