



786 Swing Head

2.786.0040

納品範囲 2.786.0040

Qt.	Order no.	説明
1 PCS	1.786.0010	786 Swing Head with transfer head left



1 PCS

6.2621.120

Hexagon key 1.5 mm



オプションアクセサリ

Order no.	説明	
6.1462.240	Transfer head for 786 Swing Head, bent and right swinging	
6.1462.250	Swing arm for Polytron and 786 Swing Head, left swinging	
6.1462.260	Robotic arm with holder for titration head, left swinging, external	
6.2058.040	Swing arm enhancement for Swing Head 786	

