




791 VA Detector

2.791.0020

納品範囲 2.791.0020

Qt.	Order no.	説明
1 PCS	1.791.0010	791 IC VA Detector



3 PCS

6.2120.020

Mini electrode cable / plug F - 2 mm



4 PCS

6.2704.010

Seal for VA detector cell



3 PCS

6.2704.020

Seal for VA detector cell



1 PCS

6.2813.020

Dummy Cell for VA Detector



オプションアクセサリ

Order no.	説明	
6.1204.130	Silver Electrode Tip	
6.1204.140	Gold Electrode Tip	
6.2128.170	Cable 771 - analog	
6.2802.000	Polishing set for solid-state electrodes	