



804 Ti Stand with stand

2.804.0040

納品範囲 2.804.0040

Qt.	Order no.	説明
1 PCS	1.804.0010	804 Ti Stand without stand rod



1 PCS

6.2001.060

Support stand without support rod



1 PCS

6.2013.010

Clamping ring



1 PCS

6.2016.070

Support rod / 400 mm



1 PCS

6.2021.020

Electrode holder



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.2621.130

Hexagon key 2 mm



オプションアクセサリー

Order no.	説明	
6.2001.090	Base plate on both sides	
6.2151.010	Cable MSB plug / socket	