



823 Membrane Pump Unit «aspirate»

2.823.0020

納品範囲 2.823.0020

Qt.	Order no.	説明
1 PCS	1.823.0010	823 Membrane Pump Unit



1 PCS

6.1446.160

Thermometer stopper, B-14/15, oblique
insertion



1 PCS

6.1543.170

Aspiration tip / M8 thread



1 PCS

6.1621.000

Container 10 L



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



1 PCS

6.1820.070

Coupling 4 mm / 6 mm / M8 inner



1 PCS

6.1828.000

PVDF connection nipple

