



823 Membrane Pump Unit «rinse»

2.823.0030

納品範囲 2.823.0030

Qt.	Order no.	説明
1 PCS	1.823.0010	823 Membrane Pump Unit



1 PCS

6.1621.000

Container 10 L



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



3 PCS

6.1805.060

FEP tubing / M6 / 60 cm



1 PCS

6.1808.170

Distributor for Sample Processors for rinsing and aspirating



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



1 PCS

6.1828.000

PVDF connection nipple



3 PCS

6.2740.020

Spray nozzle

