



826 pH mobile with carrying case

2.826.0110

納品範囲 2.826.0110

Qt.	Order no.	説明
1 PCS	1.826.0010	826 pH mobile



1 PCS

6.2050.000

Carrying strap for pH Meter



1 PCS

6.2716.040

Case for 826 pH Meter



1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm



1 PCS

6.1446.000

SGJ stopper / B-14/(15)



1 PCS

6.1613.020

Bottle / 25 mL / label pH 7 (empty)



1 PCS

6.1613.030

Bottle / 25 mL / label pH 4 (empty)



1 PCS

6.1614.000

Wash bottle / 250 mL



3 PCS

6.1614.030

Plastic bottle / 50 mL



1 PCS

6.2133.000

Batteries LR6, 1.5V, 4 pieces



1 PCS

6.2307.230

Buffer solutions pH 4, 7 and 9



1 PCS

6.2308.050

Electrolyte 3 mol/L KCl (50 mL)



1 PCS

6.2621.140

Hexagon key 2.5 mm



1 PCS

6.2717.000

PP beaker, 100 mL



オプションアクセサリー

Order no.	説明
6.2325.000	pHit kit


