



827 pH lab; 240 V, AUS with Primatrode

2.827.0117

納品範囲 2.827.0117

Qt.	Order no.	説明
1 PCS	6.2161.040	Plug-in power supply unit AUS 240V/6V DC



1 PCS

6.2013.010

Clamping ring



1 PCS

6.2016.050

Support rod / 300 mm



1 PCS

6.2021.020

Electrode holder



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



オプションアクセサリー

Order no.	説明
-----------	----

6.2307.230 Buffer solutions pH 4, 7 and 9



6.2325.000 pHit kit

