



833 Advanced IC Liquid Handling Dialysis Unit

2.833.0040

納品範囲 2.833.0040

Qt.	Order no.	説明
7 PCS	6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m



1 PCS

6.1831.050

PEEK capillary 0.5 mm i.d., 40 cm



1 PCS

6.1831.060

PEEK capillary 0.5 mm i.d., 1 m



1 PCS

6.2057.010

Holder for dialysis cell



1 PCS

6.2128.180

Cable 762/830 - 752/753/754/833



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.2714.010

Dialysis membrane (cellulose acetate)



1 PCS

6.2729.100

Dialysis cell



1 PCS

6.2744.000

Pressure screw PVDF



2 PCS

6.2744.010

Pressure screw 5x



2 PCS

6.2744.034

カップリング チューブニップル/UNF 10/32、2個


圧力調整ネジとポンプチューブを接続。2個。ペリスタリックポンプ付きIC装置用。



4 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device




1 PCS

6.2831.010

Plastic tweezers



オプションアクセサリー

Order no.	説明	
6.1826.060	Pump tubing (orange/yellow)	
6.2714.030	Dialysis membrane (polyamide)	
6.2755.000	Tubing cartridge	