



843 Pump Station (membrane)

2.843.0020

納品範囲 2.843.0020

Qt.	Order no.	説明
1 PCS	1.843.0020	843 Pump Station (membrane)



4 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



2 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



2 PCS

6.1828.000

PVDF connection nipple



オプションアクセサリ

Order no.	説明	
6.2141.260	Remote cable Compact Sample Changer - Dosimat plus and Titrino plus	
6.2141.300	Remote cable 843 Pump Station to Sample Processor	