




843 VA Membrane Pump Station コンパクト VAシステム

2.843.0040

843ダイヤフラムポンプステーション には、リモート信号を介して797 VAコンピュートレース によって直接制御できる内蔵式ダイヤフラムポンプが2つ装備されています。このポンプは、797 VAコンピュートレースを接続した一部自動および完全自動のVAおよびCVSシステムの測定容器の洗浄および排出を自動で行います。電圧電流法のバージョンでは、測定容器の自動洗浄および自動排出のための付属品が完全装備されています。

納品範囲 2.843.0040

Qt.	Order no.	説明
1 PCS	1.843.0020	843 Pump Station (membrane)



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1456.210

Measuring vessel / 10-90 mL



1 PCS

6.1602.115

Bottle attachment / S 40 / green



2 PCS

6.1621.000

Container 10 L



1 PCS 6.1805.020 FEP tubing / M6 / 52 cm



1 PCS 6.1805.100 FEP tubing / M6 / 40 cm



1 PCS 6.1805.180 FEP tubing / M6 / 16 cm



4 PCS 6.1805.530 FEP tubing / M6 / 2 m



1 PCS

6.1808.000

Connection bushing M6



2 PCS

6.1819.010

PTFE tube



4 PCS

6.1820.010

Adapter for 1/8" gas line/M6 thread



1 PCS

6.1828.020

Five-way M6 tubing connector for canister



1 PCS

6.1829.020

FEP aspiration tube M6, 0.5 m



1 PCS

6.2141.280

Connection cable 843 Pump Station to 797 VA
Computrace



1 PCS

6.2621.000

Adjustable wrench



オプションアクセサリ

Order no.	説明	
6.1621.100	Cap for 10-L container	
6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m	
6.1808.010	Coupling piece 2 X M6	
6.1812.000	PTFE tubing 4 mm / 6 mm / 4 m	

6.1820.020

M10/M6 adapter



6.1828.000

PVDF connection nipple

