



# 843 Pump Station (membrane) - rinse /aspirate for Sample Processors

2.843.0050

納品範囲 2.843.0050

| Qt.   | Order no.  | 説明                          |
|-------|------------|-----------------------------|
| 1 PCS | 1.843.0020 | 843 Pump Station (membrane) |



1 PCS

6.1543.170

Aspiration tip / M8 thread



---

2 PCS

6.1621.000

Container 10 L



---

2 PCS

6.1621.100

Cap for 10-L container



---

4 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



3 PCS

6.1805.060

FEP tubing / M6 / 60 cm



---

3 PCS

6.1805.420

FEP tubing / M6 / 48 cm



---

1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm



---

1 PCS

6.1808.170

Distributor for Sample Processors for rinsing and aspirating



2 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



---

2 PCS

6.1828.000

PVDF connection nipple



---

1 PCS

6.2141.300

Remote cable 843 Pump Station to Sample Processor



---

3 PCS

6.2740.020

Spray nozzle



オプションアクセサリー

| Order no.  | 説明             |
|------------|----------------|
| 6.2740.030 | Rinsing nozzle |



---