



843 Pump Station (peristaltic) - rinse /aspirate for Compact Sample Changer

2.843.0130

納品範囲 2.843.0130

Qt.	Order no.	説明
1 PCS	1.843.0120	843 Pump Station (peristaltic)



1 PCS

6.1543.120

Aspiration tip / M8 thread



2 PCS

6.1621.000

Container 10 L



2 PCS

6.1621.100

Cap for 10-L container



4 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1805.200

PTFE tubing / M8 / 0.5 m



3 PCS

6.1805.420

FEP tubing / M6 / 48 cm



1 PCS

6.1808.240

Distributor Rinsing / Aspiration for 869
Compact Sample Changer



2 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



4 PCS

6.1820.050

M6 / Olive adapter



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m



2 PCS

6.1828.000

PVDF connection nipple



3 PCS

6.2740.020

Spray nozzle



オプションアクセサリ

Order no.	説明	
6.2141.260	Remote cable Compact Sample Changer - Dosimat plus and Titrino plus	
6.2141.300	Remote cable 843 Pump Station to Sample Processor	
6.2740.030	Rinsing nozzle	
6.2761.000	Splash protection for 843 Pump Station (Peristaltic)	