



843 VA Peristaltic Pump Station for Compact VA systems

2.843.0140

納品範囲 2.843.0140

Qt.	Order no.	説明
1 PCS	1.843.0120	843 Pump Station (peristaltic)



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1456.210

Measuring vessel / 10-90 mL



1 PCS

6.1602.115

Bottle attachment / S 40 / green



2 PCS

6.1621.000

Container 10 L



1 PCS 6.1805.020 FEP tubing / M6 / 52 cm



1 PCS 6.1805.100 FEP tubing / M6 / 40 cm



1 PCS 6.1805.180 FEP tubing / M6 / 16 cm



4 PCS 6.1805.530 FEP tubing / M6 / 2 m



1 PCS

6.1808.000

Connection bushing M6



2 PCS

6.1819.010

PTFE tube



4 PCS

6.1820.020

M10/M6 adapter



4 PCS

6.1820.050

M6 / Olive adapter



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m



1 PCS

6.1828.020

Five-way M6 tubing connector for canister



1 PCS

6.1829.020

FEP aspiration tube M6, 0.5 m



1 PCS

6.2141.280

Connection cable 843 Pump Station to 797 VA Computrace



オプションアクセサリ

Order no.	説明	
6.1621.100	Cap for 10-L container	
6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m	
6.1808.010	Coupling piece 2 X M6	
6.1812.000	PTFE tubing 4 mm / 6 mm / 4 m	

6.1828.000

PVDF connection nipple

