



843 Pump Station (peristaltic) - rinse /aspirate for Sample Processors

2.843.0150

納品範囲 2.843.0150

Qt.	Order no.	説明
1 PCS	1.843.0120	843 Pump Station (peristaltic)



1 PCS

6.1543.170

Aspiration tip / M8 thread



2 PCS

6.1621.000

Container 10 L



2 PCS

6.1621.100

Cap for 10-L container



4 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



3 PCS

6.1805.060

FEP tubing / M6 / 60 cm



3 PCS

6.1805.420

FEP tubing / M6 / 48 cm



1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm



1 PCS

6.1808.170

Distributor for Sample Processors for rinsing and aspirating



2 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



4 PCS

6.1820.050

M6 / Olive adapter



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m



2 PCS

6.1828.000

PVDF connection nipple



1 PCS

6.2141.300

Remote cable 843 Pump Station to Sample
Processor



3 PCS

6.2740.020

Spray nozzle



オプションアクセサリー

Order no.	説明
6.2740.030	Rinsing nozzle


