



Professional VA/CVS システム用 843 Peristaltic Pump Station

2.843.0230

843 Peristaltic Pump Station は、半自動および完全自動の Professional VA/CVS システムの測定容器の洗浄および排出を自動で行います。これらには、リモート信号によって制御できる内蔵式ペリスタリックポンプが2つ装備されています。

VA と CVS のバリエーションには、測定容器の自動洗浄と排出用の付属品一式があります。接続ケーブルは別途ご注文ください。

納品範囲 2.843.0230

Qt.	Order no.	説明
-----	-----------	----

1 PCS

1.843.0120

843 Pump Station (peristaltic)



1 PCS

6.1446.000

SGJ stopper / B-14/(15)



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1602.115

Bottle attachment / S 40 / green



2 PCS

6.1621.000

Container 10 L



2 PCS

6.1621.100

Cap for 10-L container



4 PCS

6.1805.530

FEP tubing / M6 / 2 m



1 PCS

6.1819.110

FEP aspiration tubing to canister



4 PCS

6.1820.020

M10/M6 adapter



4 PCS

6.1820.050

M6 / Olive adapter



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m



1 PCS

6.1828.020

Five-way M6 tubing connector for canister



1 PCS

6.1829.070

Tubing set for Professional VA/CSV instruments



1 PCS

6.2053.000

Cable clip



オプションアクセサリ

Order no.	説明
6.2148.010	Remote Box MSB



6.2141.280	Connection cable 843 Pump Station to 797 VA Computrace
------------	--



6.2141.300	Remote cable 843 Pump Station to Sample Processor
------------	---

