



845 Eluent Synthesizer with two Dosinos

2.845.0020

納品範囲 2.845.0020


Qt.	Order no.	説明
1 PCS	1.845.0020	845 Eluent Synthesizer with two Dosinos



1 PCS

6.1602.180

Bottle Adapter to 845



2 PCS

6.1838.000

FEP Tubing Connection to 845, 1 x M6 / 45 cm



2 PCS

6.2248.010

Magnetic shield , 15 x 68 mm, 3 pieces



1 PCS

6.1546.030

エキストラクターロッド

807 Dosing Units または シリンダーユニット OMNIS のピストン用。



1 PCS

6.1608.070


Eluent bottle / 2 L / GL 45



5 PCS

6.1803.130

Capillary 1/8 in. / 1/16 in. i.d. / 1.5 m



1 PCS

6.1816.020

Silicone tubing 6 mm i.d., 1 m



1 PCS

6.1903.060

Stirring bar / 30 mm



1 PCS

6.2151.000

Cable USB A – mini-DIN 8-pin



1 PCS

6.2739.000

Wrench



2 PCS

6.2744.170

Pressure screw PEEK UNF 1/4-28, 3 pieces



1 PCS

6.2803.010

Grease (2 g)



オプションアクセサリ

Order no.	説明	
6.5329.000	Equipment for pH Adjustment at the 845 Eluent Synthesizer	
6.1602.160	溶離液および試薬のためのボトルキャップGL 45 溶離液ボトル、試薬ボトルおよび洗浄ボトル用、乾燥管および吸引チューブ用コネクタ付き。	
6.1608.030	Round glass bottle / 1000 mL / GL 45	
6.1609.000	Adsorbing tube, large and bent	

