



871 Advanced Bioscan

2.871.0010

納品範囲 2.871.0010

Qt.	Order no.	説明
1 PCS	1.871.0010	871 Advanced Bioscan



1 PCS

6.2156.010

Recorder cable for Bioscan



1 PCS

6.1816.060


Drain tubing for UV/VIS Compact IC and Bioscan



1 PCS

6.1825.210

PEEK sample loop 20 µL



1 PCS

6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m



1 PCS

6.2620.150

Pulsation absorber (metal-free)



2 PCS

6.2744.010

Pressure screw 5x



1 PCS

6.2744.014

Pressure screw 2x



1 PCS

6.2744.040

2 x UNF 10/32 coupling



オプションアクセサリ

Order no.	説明
6.1254.070	Silver flow cell for Bioscan




6.1254.130	Platinum Variocell for Bioscan
------------	--------------------------------



6.1254.140	Silver Variocell for Bioscan
------------	------------------------------



6.1254.230	Platinum working electrode for Variocell
------------	--



6.2744.130

Pressure screw for Bioscan measuring cell



6.1803.030

PTFE capillary 0.5 mm i.d. / 3 m

