



943 Professional Reactor Vario

2.943.0110

943 プロフェッショナル Reactor Varioは、特にプレカラム誘導体化およびポストカラム誘導体化のために開発された加熱できるリアクターです。この製品は、頑丈さ、高速度加熱、150℃までの反応を実行できることを特長としています。

納品範囲 2.943.0110

Qt.	Order no.	説明
1 PCS	1.943.0010	943 Professional Thermostat / Reactor Vario



1 PCS 6.1015.100 Metrosep BP 1 Guard/2.0

Metrosep BP 1 Guard/2.0は、ポストカラム試薬の流路に充分な高い作業圧力を生成するのに使用されます。これは搬送用高圧ポンプと組み合わせることで、専門リアクターでのポストカラム試薬の供給に使用されます。



1 PCS 6.1803.150 PTFE Capillary i.d. 0.25 mm / 1 m



2 PCS 6.1831.100 PEEK capillary 0.25 mm i.d., 1 m



1 PCS 6.2151.020 Cable USB A - USB B 1.8 m



1 PCS

6.2621.100

3 mm hex key for IC Sample Processors



1 PCS

6.2621.140

Hexagon key 2.5 mm



1 PCS

6.2744.010

Pressure screw 5x



1 PCS


6.2744.014


Pressure screw 2x






オプションアクセサリ

Order no.	説明	
6.1602.200	Bottle cap 2 x M6, 1 x M8 to 6.1608.120	
6.1608.120	PP eluent bottle, 2 L	
6.1836.020	Preheating capillary with core, 1.44 m	
6.2061.110	Tray with sensor for Professional IC Instruments	



6.2845.600 Column holder complete to Professional Thermostat



6.2061.100 ボトルホルダー Professional IC 装置

Professional IC 装置のための装置カバーおよびボトルホルダー。



6.2845.210 HT reactor, complete, for 6.2845.110

