



Autolab IMP

3500003360


納品範囲 3500003360

Qt.	Order no.	説明
1 PCS	CABLE.PWR	Power cable




1 PCS

CELLCBL.M101. Cell cable
204



1 PCS

AUT.DUMCELL. Autolab dummy cell
S




1 PCS

CBL.USB Standard USB cable



オプションアクセサリ

Order no.	説明
CBL.MONIT. MAC.204.S	Monitor cable for M101/M204/PGSTAT204



ALL.CLIP.BLACK Black Alligator Clamp



ALL.CLIP.RED Red Alligator Clamp



AUT204.CASE Field case for PGSTAT204



