



マイクロ電極 (長さ 17.8 cm)

6.0234.110

水性酸/塩基滴定のための細い複合pH電極、直径 6 mm、取り付け長さ 16.8 cm。
参照内部液: $c(\text{KCl}) = 3 \text{ mol/L}$ 、保存液で保管。

オプションアクセサリ

Order no.	説明
6.2308.020	Electrolyte 3 mol/L KCl (250 mL)



6.2323.000 Storage solution



6.2325.000 pHit kit



6.2325.100 洗浄液 3 x 50 mL



6.2104.020 Electrode cable / 1 m / F

