



Electrode shaft

6.1241.060

納品範囲 6.1241.060

Qt.	Order no.	説明
1 PCS	6.1236.020	Sleeve with SGJ 14/12 mm



オプションアクセサリ

Order no.	説明	
6.1204.130	Silver Electrode Tip	
6.1204.140	Gold Electrode Tip	
6.1204.150	Lateral Gold Electrode Tip	
6.1204.170	3 mm platinum electrode tip for CVS	

6.1204.180 Ultra Trace electrode tip



6.1204.190 1 mm platinum electrode tip for CVS



6.1204.610 2 mm platinum electrode tip for CVS



6.2103.110 Contact pin M3, 2 mm



6.2103.120

Contact pin M4 / 2 mm



6.1204.600

Glassy carbon electrode tip, shaft made of glass.

